



**HOSPITAL BEDS &
HOSPITAL FURNITURE**

GRACE MEDICAL INSTRUMENT CO.,LTD

TEL:+86 193 3372 3988

E-mail:supplier@gracemedy.com

Website:www.gracemedy.com

Address:No. 10 Xiwaihuan Road, Zaoqiang County
Hengshui City, Hebei Province, China



Official website



WhatsApp

**MEDICAL DEVICE MANUFACTURE
ONE-STOP SERVICE**

GRACE MEDY **PRODUCT BROCHURES**

Hospital Beds & Hospital Furniture

GRACE MEDY

Embrace The Future

Gracemedy Is Here

Escort Your Health and Spread Love



CONTENT

01	GRACE STRENGTH
03	MEDICAL BED
54	EXAMINATION BED

59	MATTRESS
62	BEDSIDE CABINET
65	DINING TABLE

67	SCREEN
69	OPERATING TABLE
75	LAMP

82	PATIENT TROLLEY
89	STRETCHER
93	MEDICAL TROLLEY
103	CHAIR



GRACE STRENGTH



Grace Medy is a professional manufacturer and designer of hospital beds and accessories, wheelchair and relate medical products. Our factory established in 2009, after years of development, we have created a professional R&D team focusing on electric ICU beds, home care beds, electric wheelchair and all the related accessories. Our comprehensive and intelligent production lines provide you with modern medical and nursing combined products.

Our company has passed CE certification, FDA, and ISO13485 quality management system certification. These certificates are a recognition of our quality.

Our commitment: First-class products, First-class quality, First-class service





BCEM-JTS05-1

 Multifunctional Hospital Electric ICU Bed

Bed length	2110-2310mm
Bed width	980mm
Bed height	620-870mm
Bed board length	1910-2110mm
Bed board width	820mm



BCEM-JTS05-1

 Multifunctional Hospital Electric ICU Bed



BCEM-JTS05-2

 Multifunctional Hospital Electric ICU Bed



BCEM-JTS05-3

 Multifunctional Hospital Electric ICU Bed

Accessories

- Stainless steel four-hook infusion pole
- 4 carbon steel drainage hooks
- Controller grounding wire
- With battery

■ Configuration

- European style removable PP head and foot
- PP four-piece European style guardrail with angle display and guardrail controller
- Back anti-bate bed board other plastic steel bed board
- Low-carbon square tube bed frame
- 6-inch central control chassis (side brake)
- Back board with CPR
- 4 Φ75 anti-collision wheels
- With middle directional wheel
- First aid push handrail
- LINAK motor



FUNCTION

- Back section adjustable: $70^{\circ} \pm 2^{\circ}$
- Knee section adjustable: $20^{\circ} \pm 2^{\circ}$
- Overall lift: 620-870mm (deviation ± 10 mm)
- Forward tilt angle: $18^{\circ} \pm 2^{\circ}$
- Backward tilt angle: $18^{\circ} \pm 2^{\circ}$
- Left and right tilt: $0-15^{\circ} \pm 3^{\circ}$
- Leg extension bed board folding down: $65^{\circ} \pm 2^{\circ}$
- Leg extension bed board can be extended 200mm (deviation ± 10 mm)
- One-button sitting position
- One-button flattening
- One-button CPR
- Backboard retraction abdominal decompression function

BCE5-RK02-1

Five-function electric ICU hospital bed

Specification

Motorized back raise, Knee raise, Height adjustment, (Reverse) Trendelenburg, Central locking, One-button sitting position, One-button flattening, One-button CPR, Backboard retraction abdominal decompression function



Color 1



Color 2



Color 3



Color 4



Color 5



Color 6



Color 7



Color 8

Standard accessory

- Electric Motors with Backup battery (Power supply: 110-220 VAC)-----4pcs
- CPR cardiopulmonary resuscitation position with Handle-----1pc
- CPR control with removable headboard and footboard-----2pcs
- Center controlled lock-----4pcs
- Foldaway railing made of ABS Plastic-----4pcs
- Urine hooks-----4pcs
- I.V. drip holes-----4pcs
- I.V. pole-----1pc

Overall size: 2150*1050*460-720mm

Bed board: 1950*830mm

Back-rest lifting angle is 0-75°(±2%);

Knee-rest lifting angle is 0-45° (±2%);

Height adjustable from 460 to 720mm;

Trendelenburg: 0~12°; Anti-Trendelenburg: 0~12°

BCE5-SMJ02-5

 multifunction weighing ICU hospital bed

Specification

multifunction hospital bed have backrest, leg rest, height adjustment, trendelenburg, weighing function and reverse trendelenburg adjustment functions.



Overall size: 2200*1040*(420-680)mm

Back section adjustable: 0-75°

Knee section adjustable: 0-45°

Adjustable height: 420~680mm

Trendelenburg / Reverse Trendelenburg: 0-15°

Color



Configuration



Headboard/Footboard
Detachable head & foot board.



Brake system
Central brake central control casters



Damping lifting side rail
Easy quick open by one touch lever



Bed surface
The parts are made by steel and environmentally materials

Optional



Weighing system



Headboard & Guardrail



Central control and central brake casters



Infusion stand

BCE5-SMJ02-13

Five-function electric ICU hospital bed

Specification

five-function hospital bed have backrest, leg rest, height adjustment, trendelenburg and reverse trendelenburg adjustment functions.



Overall size: 2090*1040*420~680mm

Bed size: 1950*900*420~680mm

Back section adjustable: 0-75°

Knee section adjustable: 0-45°

Adjustable height: 420~680mm

Trendelenburg / Reverse Trendelenburg: 0-15°

Color



Configuration



Headboard/Footboard
Detachable head & foot board in case of emergency.



Brake system
Central control side brake



Damping lifting side rail
Easy quick open by one touch lever



Bed surface
The parts are made by steel and environmentally materials

Optional



anti-collision head



Mattress



Infusion stand

BCE5-SMJ02-3

Five-function electric ICU hospital bed

Specification

five-function hospital bed have backrest, leg rest, height adjustment, trendelenburg and reverse trendelenburg adjustment functions.



Overall size: 2200*1040*(420-680)mm

Back section adjustable: 0-75°

Knee section adjustable: 0-45°

Adjustable height: 420~680mm

Trendelenburg / Reverse Trendelenburg: 0-15°

Color



Configuration



Headboard/Footboard
Detachable head & foot board in case of emergency.



Central brake
Central control casters



Damping lifting side rail
Easy quick open by one touch lever



Bed surface
The parts are made by steel and environmentally materials


Optional



Headboard & Guardrail



Central control and central brake casters



Infusion stand



Motor brand

BCE5-SMJ02-1

Five-function electric hospital bed

Specification

five-function hospital bed have backrest, leg rest, height adjustment, trendelenburg and reverse trendelenburg adjustment functions.



Overall size: 2200*1040*(420-680) mm

Bed surface size: 1950*900mm

Back section adjustable: 0-75°

Knee section adjustable: 0-45°

Adjustable height: 420~680mm

Trendelenburg / Reverse Trendelenburg: 0-15°

Configuration

With drainage hook optional
Stainless steel double wire-removing crank
Castors: optional
Guardrails: optional

Color



Configuration



Headboard/Footboard
Detachable head & foot board in case of emergency.



Damping lifting side rail
Easy quick open by one touch lever



Luxury silent casters
Independent brake, easy to operate




Bed surface
The parts are made by steel and environmentally materials

Optional



Headboard & Guardrail



Motor brand




Central control and central brake casters

BCE5-ZR02-10

Five-function electric ICU hospital bed

Specification

Motorized back raise, Knee raise, Height adjustment, (Reverse) Trendelenburg, One-button sitting position, One-button flattening, One-button CPR Backboard retraction abdominal decompression function



BACK LIFT



LEG LIFT



BACK AND LEG LIFT TOGETHER



TRENDELENBURG



REVERSE TRENDLENBURG



OVERALL LIFT

Overall size: 2150*1050*460-720mm

Bed board: 1950*830mm

Back-rest lifting angle is 0-75°(±2%);

Knee-rest lifting angle is 0-45° (±2%);

Height adjustable from 460 to 720mm;

Trendelenburg: 0~12°; Anti-Trendelenburg: 0~12°

Standard accessory

Motor brand: optional

Castors: optional 5-inch silent casters driven into the bed frame

Headboard: hanging distance 900, optional

Guardrails: four small guardrails, optional

Color



Configuration



Headboard/Footboard

Detachable head & foot board in case of emergency.



Damping lifting side rail

Easy quick open by one touch lever



Luxury silent casters

Independent brake, easy to operate



Bed surface

The parts are made by steel and environmentally materials

Optional



Headboard & Guardrail



5-inch silent casters driven into the bed frame



Motor brand

BCE5-SMJ03-3

Five-function electric ICU hospital bed

Specification

five-function hospital bed have backrest, leg rest, height adjustment, trendelenburg and reverse trendelenburg adjustment functions.



BACK LIFT



LEG LIFT



BACK AND LEG
LIFT TOGETHER



TRENDELENBURG



REVERSE
TRENDELENBURG



OVERALL
LIFT

Overall size: 2200*960*(420-680)mm

Back section adjustable: 0-75°

Knee section adjustable: 0-45°

Adjustable height: 420~680mm

Trendelenburg / Reverse Trendelenburg: 0-15°

Color



Configuration



Headboard/Footboard

Detachable head & foot board in case of emergency.



Brake system

Central brake central control casters



6Rod aluminum alloy foldable siderail



Mesh surface

The mesh surface has better air permeability


Optional



Headboard & Guardrail



Double Side Caster With Brake



Infusion stand

BCM5-SMJ02-1

Five-function manual hospital bed

Specification

five-function hospital bed have backrest, leg rest, height adjustment, trendelenburg and reverse trendelenburg adjustment functions.



BACK LIFT



LEG LIFT



BACK AND LEG
LIFT TOGETHER



OVERALL LEAN
FORWARD



REVERSE
TRENDLENBURG



OVERALL
LIFT

Configuration

With drainage hook optional
Stainless steel double wire-removing crank
Castors: optional
Guardrails: optional

Color



Configuration



Headboard/Footboard

Detachable head & foot board in case of emergency.



Damping lifting side rail

Easy quick open by one touch lever



Luxury silent casters

Independent brake, easy to operate



Bed surface

The parts are made by steel and environmentally materials


Optional



Headboard & Guardrail



Central control and central brake casters



Infusion stand

Overall size: 2200*1040*(420-680) mm

Bed surface size: 1950*900mm

Back section adjustable: 0-75°

Knee section adjustable: 0-45°

Adjustable height: 420~680mm

Trendelenburg / Reverse Trendelenburg: 0-15°

BCM3-SMJ02-2

Three-function manual hospital bed

- Overall size: 2090*1040*420-680mm;
- Bed size: 2020*900*420~680mm
- Back section adjustable: 0-75°
- Knee section adjustable: 0-45°
- Adjustable height: 420~680mm

Standard configuration

- ✓ 5-inch casters: central brake central control casters*4
- ✓ Crank: ABS folding hidden crank
- ✓ Four small guardrails*1 set
- ✓ Head and foot boards with anti-collision wheels*1 set



Optional

- IV stand holes: 4 or 2 diagonal holes optional
- Drainage bag hooks: 2
- IV stand: optional



BCM3-SMJ02-8

Three-function manual hospital bed

- Overall size: 2090*960*420-680mm;
- Bed size: 2020*900*420~680mm
- Back section adjustable: 0-75°
- Knee section adjustable: 0-45°
- Adjustable height: 420~680mm

Standard configuration

- ✓ 5-inch casters: central brake central control casters*4
- ✓ Crank: ABS folding hidden crank (replaceable stainless steel crank)
- ✓ Guardrail (gear guardrail)*1 set
- ✓ Head and foot board*1 set



Optional

- IV stand holes: 4 or 2 diagonal holes optional
- Drainage bag hook: 2
- IV stand: optional



BCE3-ZR02-1

Three-function electric hospital bed

- Overall size: 2090*1040*420-680mm;
- Bed size: 2020*900*420~680mm
- Back section adjustable: 0-75°
- Knee section adjustable: 0-45°
- Adjustable height: 200mm

Standard configuration

- ✓ 5-inch single-sided casters*4
- ✓ Motor: 1 set of 3 motors
- ✓ Four small guardrails (with angle display)*1 set
- ✓ Head and foot boards*1 set



Optional

- IV stand holes: 4 or 2 diagonal holes optional
- Drainage bag hook: 2
- IV stand: optional



BCE3-ZR02-3

Three-function electric hospital bed

- Overall size: 2090*960*420-680mm
- Bed size: 1950*900*420~680mm
- Back section adjustable: 0-75°
- Knee section adjustable: 0-45°
- Adjustable height: 420~680mm

Standard configuration

- ✓ 5-inch single-sided casters*4
- ✓ Motor: 1 set of 3 motors
- ✓ Four small guardrails (with angle display)*1 set
- ✓ Head and foot boards*1 set



Optional

- IV stand holes: 4 or 2 diagonal holes optional
- Drainage bag hook: 2
- IV stand: optional



BCM3-ZR01-2

Three-function manual hospital bed

- Overall size: 2090*960*420-680mm
- Bed size: 1950*900*420~680mm
- Back section adjustable: 0-75°
- Knee section adjustable: 0-45°
- Adjustable height: 420~680mm

Standard configuration

- ✓ 5-inch casters*4
- ✓ Crank: stainless steel crank*3 (ABS crank can be replaced)
- ✓ Guardrail (gear guardrail)*1 set
- ✓ Head and foot board*1 set



Optional

- IV stand holes: 4 or 2 diagonal holes optional
- Drainage bag hook: 2
- IV stand: optional



BCM2-SMJ03-2

Two function manual hospital bed

- Overall size: 2090*960*420-680mm
- Bed size: 1950*900*420~680mm
- Back section adjustable: 0-75°
- Knee section adjustable: 0-45°
- Adjustable height: 420~680mm

Standard configuration

- ✓ 5-inch casters*4
- ✓ Crank: ABS folding crank*2(stainless steel crank can be replaced)
- ✓ Guardrail (gear guardrail)*1 set
- ✓ Head and foot board*1 set



Optional

- IV stand holes: 4 or 2 diagonal holes optional
- Drainage bag hook: 2
- IV stand: optional



BCE2-RT01-28

Two-function electric hospital bed

- Overall size: 2090*960*500mm
- Back section adjustable: 0-75°
- Knee section adjustable: 0-45°



Standard configuration

- ✓ 5-inch double-sided casters*4
- ✓ Motor: motor*2
- ✓ Guardrail (gear guardrail)*1 set
- ✓ Head and foot board*1 set

Optional

IV stand holes: 4 or 2 diagonal holes optional
Drainage bag hooks: 2
IV stand: optional



BCM2-SMJ02-33

Two function manual hospital bed

- Overall size: 2090*960*500mm
- Back section adjustable: 0-75°
- Knee section adjustable: 0-45°



Standard configuration

- ✓ 5-inch casters with side brake system
- ✓ Crank: ABS folding hidden crank*2(Stainless steel crank can be replaced)
- ✓ Guardrail (gear guardrail)*1 set
- ✓ Head and foot board*1 set

Optional

IV stand holes: 4 or 2 diagonal holes optional
Drainage bag hook: 2
IV stand: optional



BCM2-SMJ02-31

Two function manual hospital bed



- Overall size: 2200*1040*550mm;
- Back section adjustable: 0-75°
- Knee section adjustable: 0-45°



Standard configuration

- ✓ 5-inch double-sided casters*4
- ✓ Crank: ABS folding hidden crank*2
- ✓ Four small guardrails (with angle display)*1 set
- ✓ Head and foot boards with anti-collision wheels*1 set

Optional

- IV stand holes: 4 or 2 diagonal holes optional
- Drainage bag hooks: 2
- IV stand: optional
- Stainless steel crank optional



BCM2-ZR01-2

Two function manual hospital bed

- Overall size: 2090*960*500mm
- Bed size: 1950*900*500mm
- Back section adjustable: 0-75°
- Knee section adjustable: 0-45°



BCM1-WS01-22

One-function manual hospital bed

- Overall size: 2090*960*500mm
- Back section adjustable: 0-75°



BCM1-SMJ02-5

One-function manual hospital bed

- Overall size: 2090*960*500mm
- Back section adjustable: 0-75°

Standard configuration

- ☒ Guardrail *1 set
- ☒ Head and foot board*1 set



BCM1-SMJ01-9

One-function manual hospital bed

ETM2-DF01-1

Two-function manual children hospital bed



Specifications

Overall size: 1890*900*1260mm
Back section adjustable: 0-75°
Knee section adjustable: 0-45°

Accessories

Casters*4pcs
Fully surrounded guardrail*1set
Head and footboard*1set



BCEM-MY05-1

One-function manual children hospital bed



Overall size: 1880*860*500mm
Back section adjustable: 0-75° (±5°)

Accessories

- ✓ Guardrail*1 set
- ✓ Head and foot board*1 set
- ✓ Crank: ABS folding hidden crank
- ✓ Optional: casters



BCM1-SMJ-03

One function children's hospital bed

- Overall size: 1880*860*500mm
- Back section adjustable:0-75° (±5°)

- Accessories
- ✓ Guardrail*1 set
 - ✓ Head and foot board*1 set
 - ✓ Crank: ABS folding hidden crank
 - ✓ Optional: casters



YEME-PSDAB-4

Baby cot

- Overall size: 840mm+500mm780-980mm
- Back section adjustable:0-12°



YEM1-SDPS-1

Baby cot

- Overall size: 840*500*780-980mm
- Back section adjustable:0-12° (±5°)



BCEM-SMJ03-1

Wheelchair bed

Specification

- Overall size: 2180*1220*490mm
- Bed size: 2000*1000mm
- Wheelchair size: W590mm, H490mm
- Handrail to ground height 1020mm
- Pocket hole size: 320*210mm



BACK LIFT



LEG LIFT



BACK AND LEG LIFT TOGETHER

Standard accessory

- Double-sided nylon casters for bed, silent wheels for wheelchair
- Motor: motor*1 set
- Guardrails (two large guardrails)*1 set
- Head and foot boards*1 set
- IV stand holes: 4 or 2 diagonal holes optional
- 8cm sponge mattress*1
- Optional: IV stand

GRACE MEDY

BCMM-ZR05-1

Manual nursing bed

Accessories	
☒	5-inch casters*4
☒	Motor: motor*1 set (with controller*1, handle*1)
☒	Guardrail (gear guardrail)*1 set
☒	Head and foot board*1 set
☒	IV stand holes: 4 or 2 diagonal holes optional
☒	Drainage bag hook: 2
☒	IV stand: *1
☒	Stand up puller*1 set
☒	Bed pan*1
☒	Wash basin*1
☒	Dining table*1
☒	8cm mattress*1

- Overall size: 2090*960*530-695mm
- Back lift 0-75° Leg lift 0-45°



NURSING
BED

BCM3-ZR05-3

Manual full-curve nursing bed

Standard configuration	
☒	· 5-inch casters * 4
☒	· Crank handle: stainless steel crank handle * 2
☒	· Guardrail (gear guardrail) * 1 set
☒	· Head and foot board * 1 set



- Optional**
- IV stand holes: 4 or 2 diagonal holes optional
 - Drainage bag hooks: 2
 - IV stand: optional
 - Bed pan: optional
 - Shampoo basin: optional
 - Mattress: optional

- Overall size: 2050*960*560mm
- Back lift 0-75° (±5°)
- Turning over 0-45° (±5°)
- Leg lift +25-80° (±5°)



NURSIN
BED

BCZL-YLL-1

Standing rehabilitation training bed

Accessories	
Two ABS guardrails*1 set	Overall size: 2130*1260*540±20mm
Head and foot boards*1 set	Bed size: 1960*900±20mm
Foot support board*1 piece	Foot pedal size: 490*300±20mm
Training table*1 piece	With emergency stop function
Fixed belt*4 pieces	One-button reset function
Casters*4 pieces	
Motor*1 set (with hand controller)	



TRAINING
BED



BCZL-YLL-3

Standing rehabilitation training bed

Overall size: 1900*800*575±20mm
Bed size: 1800*545±20mm
Footboard size: 520*300±20mm
Standing function: Angle 0-80°±5°
With emergency stop function

Accessories
Footrest*1
Training table*1
Fixed belt*4
Upper limb trainer*1
Casters*4
Motor*1
Controller*1



TRAINING
BED



BCZL-YLL-2

Standing rehabilitation training bed

- Overall size: 1900*800*575±20mm
- Bed size: 1775*520±20mm
- Foot pedal size: 490*300±20mm
- Manual back lifting angle: 0-65°±10°
- Standing angle: 0-80°±5°
- With emergency stop function

Accessories	
✓	Guardrail*1 set
✓	Foot support board*1 piece
✓	Training table*1 piece
✓	Fixed belt*4 pieces
✓	Upper limb trainer*1 piece
✓	Casters*4 pieces
✓	Motor*1 set
✓	Manual rocker*1 piece



NURSING
BED

BCEM-DB05-1

Intelligent electric patient transfer trolley



- Overall size: 2020*710*600-870mm
- Bed length: 180mm
- Bed width: 650mm
- Board width: 450mm (allowable deviation ±20mm)
- Bed side thickness: 1.2mm
- Bed side square tube thickness: 40×60mm; 30×60mm
- Main bed connection parts: 3mm cold plate stamping
- Bed surface: thickened high-quality stainless steel plate
- Bed surface load-bearing weight: ≥200kg

Supporting items


- controller, battery, motor, whole machine frame, manual, warranty card



HLE5-XLL-1

Electric five-function nursing bed

- Overall size: 2170*1000* (335-725) mm
- Back lifting angle: 0~70±5°
- Leg lifting angle: 0~20±5°
- Head and tail tilt angle: 12±2°



NURSING
BED

HLE5-RQ05-1

Electric five-function nursing bed

- Overall size: 2200*1030*330-730 mm
- Back lifting angle: 0~70±5°
- Leg lifting angle: 0~20±5°
- Head and tail tilt angle: 12±2°



A.02.02.03.12

Electric nursing bed


- Overall size: 2120x960x500mm
- Bed board: Cold-roll steel,podwer coating
- Head and rear bed board: Plastic
- Guardrail: Auminum alloy,collapsible type
- Castor: Universal caster, two casters with brake
- With a bed pan,lift pole



HLQUE-ZR-3

Electric nursing bed

- Overall size: 2090*960*530~695mm
- Anti-slip function: Prevent mattress from slipping and sliding
- Small wheel assistance: full-bend floor stand with small wheels



HLE3-CJ03-2

Electric three-function nursing bed

- Overall size: 2080*980*420-700mm
- Height adjust: 420-700mm
- Back lifting angle: 0-70°±5°
- Leg lifting angle: 0-30°±5°



HLE3-CJ03-3

Electric three-function nursing bed

- Overall size: 2080*1040*450mm
- Back folding angle ≥70°,
- Leg folding angle ≥30°,
- Leg bending angle ≥50°



HLM2-CJ05-1

Two-crank manual nursing bed

- Overall size: 2060*980mm
- Height from bed surface to ground: 450mm
- Back lifting angle: 0-70° ±5°
- Leg lifting angle: 0-30° ±5°



ETM2-SMJ01-1

Children orthopedic bed

- Overall size: 1880*860*500mm
- Rear lifting angle: 0-75° (±5°)
- Leg lifting angle: 0-45° (±5°)

Accessories	
☒	Handle: ABS handle*2 pcs
☒	Guardrail (gear guardrail)*1 set
☒	Head and foot board*1 set
☒	Orthopedic traction frame with lifting handle*1 set
☒	Optional: casters



ORTHOPEDIC
BED

BCMGK-SMJ01-01

Orthopedic Bed

Specification

- Overall size: 2090*960*500mm
- Bed surface length: 1950mm
- Bed surface width: 900mm
- Rear lifting angle: 0-75° (±5°)
- Leg lifting angle: 0-45° (±5°)

Accessories	
☒	Castors*4
☒	Handle: Stainless steel handle*2
☒	Guardrail (gear guardrail)*1 set
☒	Head and foot board*1 set
☒	Stainless steel gantry traction frame*1 set
☒	Optional: low pressure rod, drain, shelf



BCMGK-SMJ01-02

Orthopedic Bed

- Overall size: 2030*960*500mm
- ID length: 1950mm
- ID width: 900mm
- Back lifting angle: 0-75° (±5°);
- Leg lifting angle: 0-45° (±5°)



JCE3-KYPS-1

Electric three-function examination bed



- Overall size: 1910*600*700-950mm
- Adjustable height: 700-950mm
- Back tilt: 60±5°
- Leg tilt: 60±5° up; 30±5° down

- Standard configuration
- | | |
|--|---|
| ☒ 1 high-grade medical motor | ☒ 1 control box |
| ☒ 1 foot controller | ☒ 4-inch luxury casters |
| ☒ 4 floor brakes | ☒ 1 6cm thick mattress |



JCM1-KXDSS-1

Gynecological examination bed

- Overall size: 1850*600*820mm
- Material: #201 stainless steel
- Thickness: 1.0mm
- Mattress thickness: 60mm (including plywood)
- Mattress color: pink, black
- Function: back, hips, and legs are adjustable
- Style: assembled



- | Accessories | |
|-------------------------------------|--|
| ☒ | Freely sliding armrests * 1 pair |
| ☒ | Ordinary footrests (fixed angle, adjustable height) 1 pair |
| ☒ | Dirty basin * 1 (square, round optional) |



FKM2-BXSS-1

Gynecological examination bed

- Overall size: 1800*600*750mm
- Color: black
- Accessories: shoulder rest, leg rest



FKM2-BKPS-1

Gynecological examination bed

- Overall size: 1680*600*800mm
- Back angle adjustable
- Leg angle adjustable
- Color: blue/pink optional



JCM1-BKPS-3

Single-crank examination bed

- Overall size: 1900*600*650mm
- Bed frame: 40×40mm rectangular square tube, thickness 0.8mm
- Bed legs: 40×40mm rectangular square tube, thickness 0.8mm
- Back angle: ≥75°
- Leather color: blue
- Sponge thickness: 6cm



EXAMINATION BED

JCM0-BKPS-02

Examination bed

- Overall size: 1800×600mm×650mm
- Material: Steel spray-coated
- Bed frame 20*40mm raw material, 0.8mm thickness;
- Bed legs 40*40mm, 0.8mm thickness
- Mattress: Mattress thickness can be customized
- Color: Light brown



JCM0-BKPS-05

Examination bed

- Overall size: 1870*600*700mm
- Bed surface: PU surface, filled with 4cm sponge, 1cm wooden board
- Bed legs: 201 stainless steel, wall thickness 0.8mm
- Weight: 20kg



JCM1-BKPS-02

Examination bed


- Overall size: 1800*600*650mm
- Color: pink
- Handle: 1 stainless steel handle (or replaceable pull bar)
- Headrest *1



JCM1-BKSS-2

Examination bed

- Overall size: 1900*600*650mm
- Iron pipe thickness: 0.7mm
- Frame 25×38mm, legs 38*38mm
- Lifting method: chain lifting
- Total height: 65cm
- Mattress: 4cm sponge black mattress
- Backboard should be added with black cloth



CDDG-SMJ-3

Fully flip smart care mattress

- Overall size: 1920*960*100mm
- Function: Smart nursing mattress, no need to change the bed, one-button back lift, leg bend and turn over
- Plug-in guardrail*1 set
- Mattress*1 piece
- Motor*1 set
- Dining table: optional



GCY-FD-1

Patient slide transfer board

- Overall size: 1730*500*18mm
- Folding size: 850*500*36mm



CDQC-JHQ-1



CDQC-JHQ-2



CDQC-CTT-7





CDQC-JHQ-1

Anti-bedsores air mattress

CDCS-FLM08-1

Medical mattress


- size: can be customized according to needs
- Material: Oxford cloth, PU, composite PU optional
- Thickness: can be customized according to needs
- Lining: brown silk, sponge, sponge + palm silk



CDCS-JY08-1

Mattress

- Size:1930*900mm
- Thickness:6 cm, 8 cm, 10 cm optional
- Fabric: PU waterproof fabric
- Padding: palm fiber+sponge, polyester, high density sponge, high resilience memory foam
- Function: 1 function, 2 functions, toilet hole, optional
- Customization: Thickness, Fabric, Size, Filling



CDCS-FLM08-1

ABS/PP bedside cabinet



ModelNo:G-1G

- Overall size: 480*480*760mm
- Material: ABS or PP
- Remarks: Stainless steel finish, casters, thermos rack, shoe rack are optional

ModelNo:G-2A / G-2B(stainless steel veneer)

- Overall size: 450*420*740mm
- Material: ABS or PP
- Remarks: packing with parts or end product

ModelNo:G-2C(casters)

- Overall size:450*420*780mm
- Material: ABS or PP
- Remarks: packing with parts or end product



HLQUE-ZR-2

ABS/PP bedside cabinet

- Overall size: L490*480*780mm
- Material:ABS or PP
- Remarks: The 5 faces of the tabale are molded in one time.no joint close,steel frame inside,high strength.



CTAG-WKM-3

ABS carbon steel bedside cabinet

- Overall size: 465*455*810mm
- Material: ABS and carbon steel



G-4A



G-4B

BCMM-ZR05-1

ABS/PP bedside cabinet

Model No:G-4A

- Overall size: 480*480*780mm
- Material: ABS

Model No:G-4B

- Overall size:480*480*780mm
- Material: ABS
- Remarks: Two drawers,the food and themedicine are placed separately



CTAG-WKM-5

ABS/PP bedside cabinet

- Overall size: 502*502*765mm
- Material:ABS




CZP-LYCB-1

Extendable dining table

Specifications:

- Size: 880-1060*320*25mm
- Material: ABS
- Features: retractable, waterproof, with hook



CZAB-WKYD-1

ABS carbon steel bedside cabinet

Specifications:

- Desktop size: 820*450mm
- Height: 840-1080mm
- Main body material: carbon steel
- Load bearing: 30kg



CZAB-YKCB-2

Extendable dining table

Specifications:

- Size: 880-1060*320*25mm
- Material: ABS
- Features: retractable, waterproof, with hook



CZMZ-JLYD-1

Liftable dining table

Specifications:

- Desktop size: 765*380mm
- Height 680-1000mm
- Main body material: carbon steel
- Load bearing: 20kg



CZMZ-WKYD-1

Wooden mobile dining table

Specifications:

- Desktop size: 900*400mm
- Height: 750-1140mm
- Main body material: carbon steel
- Load bearing: 30kg



CZMZ-JLYD-2

Liftable dining table

Specifications:

- Desktop size: 600*400mm
- Height 650-1000mm
- Main body material: carbon steel
- Load bearing: 20kg



PFSS-BK4-2

Four-fold medical partition screen

Specification

- Material: Stainless steel
- Features: Removable and washable
- Height 180*Width 50cm (width/folded)
- Castors: Iron wheels with brakes
- Thickness: 0.5mm



SSE3-SXSS-1 Operating table

Specification

- Tabletop size: 2100 (±20)*500 (±20)*750 (±20) mm
- Height adjustment 250 mm
- Tilt ≥ 18°
- Forward tilt ≥ 30°
- Reverse forward tilt ≥ 15°
- Backboard adjustment ≥ 70°-18°
- Headboard adjustment ≥ ± 45°-90°
- Kidney bridge adjustment ≥ 120 mm
- Leg board ≥ 90°



SSE6-SDYKL-1 Hydraulic operating table

Specification

- Headboard: fold up ≥40°, fold down ≥90°, detachable, manual
- Backboard: fold up ≥75°, fold down ≥0°, electric
- Waist bridge: lift ≥80mm, manual
- Forward and backward tilt: ≥20°, ≥15°, electric
- Left and right tilt: ≥20° ≥15°, electric
- Longitudinal translation: 300mm±20mm, electric
- Legboard: abduction ≥90°, detachable, manual
- Tabletop: 2000mm×500mm±20mm
- Tabletop lifting: 730mm--980mm±20mm, electric
- Motor: with battery



SSE5-SDYKL-1 Hydraulic operating table

Specification

- Headboard: Up fold ≥ 45 ° Down fold ≥ 90 ° Removable manual
- Backboard: Up fold ≤ 75 ° Down fold ≥ 10 ° Electric
- Waist bridge: lift ≥ 100mm manually
- Forward and backward tilt ≥ 28 ° ≥ 18 ° electric
- Left and right tilt ≥ 23 ° ≥ 18 ° electric
- Vertical translation 400mm ± 20mm electric
- Legboard: Detachable manual with outward extension ≥ 180 °
- Tabletop: 2010mm × 500mm ± 20mm
- Table lifting: 730mm --980mm ± 20mm electric
- Power supply, power: 220V/50Hz, 250VA continuous operation with battery for 40-80 surgeries



FKEM-SDSS-2 Pulling-type electric delivery bed

Specification

- Backboard: fold up ≤60°, fold down ≥0°, electric
- Forward and backward tilt: ≥20°, electric
- Legboard: abduction ≥180°, detachable, manual
- Table size: 1200*600mm
- Auxiliary table size: 580*540mm



FKEM-SDSS-3 Gynecological maternity bed

Specification

- Backboard: fold up ≤75°, fold down ≥10°, electric
- Legboard: extension ≥180°, detachable, manual
- Tabletop: 1850*600±20mm



FKEM-SDSS-6 Gynecological delivery hospital bed



Specification

- Table size: 2020*1020mm
- Lifting height: 1000~720mm
- Table front and back tilt angle: 0-7°
- Backboard: up fold 70°
- Power supply, power: 220V/50Hz, 150VA, with battery
- Leg board: abduction = 180°

FKES-SDPS-1

Three-function electric gynecological operating table

Specification

- Backboard: fold up ≤75° fold down ≥18°
- Back tilt angle: ≥25°
- Table size:1450*590±20mm
- Height adjustment : 580--1000±20mm
- Power supply, power: 220V/50Hz, 150VA
- Power system: powered by electric push rod
- Bed frame: high-quality carbon steel spraying
- Base: leak-proof design



FKE3-SDSS-1

Three-function electric gynecological operating table

Specification

- Table size: 1900*600*(±20)mm
- Height adjustment: 650-870(±20)mm
- Forward tilt angle: ≥8°, backward tilt angle ≥20°
- Backboard adjustment angle: ≥60°-7°



YKE-QF-3

Dental Chair

Specification

- Power supply: Rated power supply: ac.220V, frequency: 50Hz, input power: 800VA
- Water heater: Rated power supply: ac.220V, frequency: 50Hz, input power: 400VA
- Water temperature constant temperature: 40°C±5°C
- Mechanical disk: Rotation angle: >270°
- Disk bearing mass ≤5KG
- Disk inclination ≤3°



Dental patient chair

- Power supply voltage:a.c. 24 V, Total load capacity: 135kg
- Air source: Air inlet pressure: 0.5MPa~0.8Mpa
- Water source: Water inlet pressure: 0.2MPa~0.4Mpa

Oral lamp

- Rated power supply: a.c.12V, power: 50VA
- Luminance: 8000 lx*20000lx

Film viewer


- Rated power supply: a.c.220V, power: 7VA
- Light source color temperature: >6500K
- Brightness: >2000cd/m²
- Stability: ≤2%
- Dispersion coefficient: ≥0.9

WYYP-SDXD-2

LED operating light with camera system



With the advancement of technology and the improvement of user requirements, only providing high-quality lighting effects can no longer meet the needs of users.



Therefore, we provide various configurations for LED series shadowless lamps, including domestic round arms, imported round arms, and imported square arms. It is also equipped with high-quality video and camera system overall solutions for users to choose from. The cameras include built-in cameras and external cameras.

Main technical parameters of LED series shadowless light camera

Camera	Lens
Sensor type : 1/2.8" Progressive Scan CMoS	Focal length : 4.7-47mm, 10x optical zoom
Scan system : line-by-line scan	Light value : F1.6-F3.0
Signal system : PAL/NSTC	Horizontal field of view: (wide-angle-telephoto)
Video output : SDI output, network output	Close-up : 120mm-800mm (wide-angle-telephoto)
Resolution : "50HZ-50fps(1920X1080)60HZ-60fps(1920x1080)50HZ-25fps(1920x1080)60HZ-30fps(1920X1080)"	Interface function
Minimum illumination : Color: 0.05Lux@(F1.6, AGC ON)	External InterfaceSDI, RS485, Power UART
Signal to noise ratioGreater than 52dB	External communicationPelco P, Pelco D, VISCA
Parameter settings	
Focus mode: semi-automatic/automatic/manual	
Exposure mode : Auto Exposure/Aperture Priority/Shutter Priority/Manual Exposure	
Electronic shutter : 1/1 second - 1/30,000 second	
White balance : Auto1/Auto2/Indoor/Outdoor/Manual/Sodium Lamp/FluoRescent Lamp	
Gain : 16 levels adjustable	
Wide dynamic : Off/On, 15 levels adjustable	
Digital noise reduction : 3D	
Lens initialization : built-in	
Image quality adjustment: Brightness, contrast, sharpness, saturation can be adjusted	
Digital zoom : 10	
General Specifications	
Power supply : DC12V±25%	
Working temperature and humidity: -10°C-60°C, less than 90%RH	
Power consumption : 5w	
Size : 43.5*45.4*68.8mm	
Weight : 190g	


WDXD-SDHB-1

Shadowless operating light

Specifications
· Illumination (Lux): 60000-160000Lux / 60000-140000Lux
· Number of bulbs: 61 / 39
· Average life of lamp beads: ≥60000h / ≥60000h
· Color temperature (adjustable): 3000K ~ 6700K / 3000K ~ 6700K
· Display index (adjustable): ≥95% / ≥95%
· Focus depth: ≥800mm / ≥800mm
· Spot diameter: ≥180mm / ≥170mm



- Very long service life
- Perfect cold-light effect
- Ideal appearance design meeting the purification requirements
- Humanized interface design
- Super deep lighting, perfect shadowless effect
- Color temperature-adjustable function



WDXD-SDHB-10

SHADOWLESS
OPERATING LIGHT

WDXD-SDYP-1

LED operating light with camera system

Specifications

- Model : The LED-700 / 500 type
- Number of LED light bulbs : 80 / 48
- Illuminance (Lux) : 60000–160000/60000–140000
- Colour temperature (K) : 3000K ~ 6700K / 3000K ~ 6700K
- Spot diameter (mm) : 150-350
- The dimming system : Inpolar dimming system
- color rendering index : ≥96
- Lighting depth (mm) : ≥1200



WYYP-SDYD-1



ZYM-MJCG-2



Patient transfer docking bed



Specification

- Product size: 365*64/73*64/9cm
- Height adjustment: 0-33cm
- Brake system: one-foot central control brake system
- Infusion stand socket: 1pc infusion stand socket at the head and foot of the bed, with a diameter of 20mm
- Drainage bag hook: equipped with four drainage hooks
- Infusion stand: stainless steel two-stage four-claw infusion stand, diameter 19mm
- Guardrail: two-piece engineering plastic damping guardrail
- Bottom cover: injection molded ABS engineering plastic bottom cover
- Back lifting angle: 0-75°
- Caster: 8pcs casters

ZYM-JHYY-1 Hydraulic transfer trolley



TRANSFER
TROLLEY

Specification

- Product size: 1930*700*630-930mm
- Back lifitng angle : 0-≥70°
- Trendelenburg / Reverse Trendelenburg: 0-≥12°
- Anti-front height: 0-≥12°
- Maximum dynamic safety load: 135kg
- Material: PP bed surface; PP guardrails;

ZYM-JHCG-1 Patient transfer trolley



Specification

- Product size: 1930*640*560-860mm
- Back lifting angle: 0-75°
- Maximum load: 250Kg
- Safe working load: 135KG
- Spray-coated square tube frame
- Bed surface: ABS bed surface
- Guardrail: ABS damping lifting large guardrail
- Caster: 150mm center control wheel
- Infusion stand: telescopic infusion stand
- Accessories: oxygen cylinder stand
- Waterproof and flame-retardant mattress: thickness 20mm
- Shockproof system: 5-level shockproof system

ZYM-MJCG-4 Conventional patient transfer trolley



Specification

- Product size: 1930*663*510-850mm
- Back angle adjustment: 0-70°
- Safe working load: 170KG
- Bed board: PP resin molded product
- Frame: Made of steel/partial aluminum products
- Casters: Central locking of 150mm diameter
- Oxygen cylinder rack: diameter 105-115mm
- Transport mattress: Waterproof fabric, easy to clean, zippers at the four corners, washable outer fabric
- Guardrail: PP resin molded guardrails on both sides
- Guardrail height: ≥290mm
- Thickness of the thinnest part of the guardrail: ≥15mm

ZYM-MJYY-4 Patient transfer trolley

Specification

- Product size: 2100*850*630/920mm
- Back lift 0 ~ 75° (±2°),
- Knee lift 0 ~ 45° (±2°),
- Tilt adjustment -15°~ 15° (±2°)
- Static load: 250KG
- Safe working load: 200KG
- Bed board: 6mm thick anti-bate material, the bed surface is a combination of three folds and four pieces
- Castors: Central control silent casters with a diameter of 200mm
- Chassis: One-piece ABS chassis, oxygen cylinders can be placed
- Chassis load-bearing: 10KG
- Guardrail: One-piece inward-retracted S-shaped guardrail, the guardrail handrail is made of aluminum alloy and 304 stainless steel



DJSS-HZXCS-1 Stainless steel patient transfer trolley



Specification

- Product size: 1800*500*700mm;
- Function: Folding leg mechanical structure, controlled by left and right handles
- Load-bearing capacity: 150kg
- Link: Stainless steel connector
- Mattress: Leather pad integrated

XZBZE-MJ-1 Hydraulic bath bed

Specification

- Overall bed size: 2150*810mm*800-1055mm
- Bed surface size 600*1950mm
- Bed surface material: 304 stainless steel plate
- Guardrail material: sus304 stainless steel
- Guardrail size: 25*25mm
- Bed body link load-bearing material: 3.0 cold-rolled steel pipe
- Hydraulic lifting system: equipped
- Front and rear tilt angle: 0~15°
- Guardrail function: Both sides of the guardrail can be rotated 80°
- Castor: 8-inch



DJAL-MDJCS-1

Aluminum alloy stretcher trolley

- Highest position 1900*550*920mm
- Lowest position 1900*550*250mm
- Maximum backrest inclination angle: 75°
- Load-bearing capacity: ≤159KG
- Connector: Aluminum alloy die-casting
- Mattress: Sponge cushion
- Material: High-quality aluminum alloy



DJAL-MDJCS-2

Aluminum alloy stretcher trolley

- Highest position: 1900*550*880mm
- Lowest position: 1900*550*320mm
- Function: with back lifting function
- Load bearing: ≤159KG
- Mattress: Oxford cloth
- Castor: Brake caster
- Material: Aluminum alloy



DJAL-MDJCS-3

Aluminum alloy stretcher trolley

- Material: Made of high-quality aluminum alloy, with aluminum alloy die-casting for the fittings
- Fixed bracket fixed to the ambulance
- Highest position: 1900*550*920mm
- Lowest position: 1900*550*250mm
- Maximum backrest inclination angle: 75°
- Load-bearing ≤159KG



DJAL-YCCS-3

Aluminum alloy stretcher trolley

- Highest position: 1930*560*840mm
- Lowest position: 1930*550*300mm
- Maximum backrest inclination: 75°
- Load-bearing capacity: 180KG
- Connectors: aluminum alloy, plastic
- Mattress: 10cm sponge cushion, PU cover
- Castors: Φ150MM125MM wide rubber wheels
- Material: aluminum alloy

DJAL-AL2-2

Aluminum alloy two-fold stretcher

- Product size: 2080*550*130mm
- Folded size: 1040*170*90mm
- Material: Aluminum alloy folding stretcher
- Wheels: with 2 wheels
- Load bearing: ≤159KG



DJDX-AL2-1

Galvanized iron two-fold stretcher

- Product size; 2000*500*150mm
- Folded size: 1070*210*110mm
- Tube material: galvanized iron
- Color: blue, camouflage, orange
- Load capacity: 120KG



DJAL-AL4-1

Aluminum alloy four-fold stretcher

- Product size: 2080*550*130mm
- Folding size: 1040*170*90mm
- Material: Aluminum alloy folding stretcher
- Weight: 4.6KG
- Load-bearing capacity: ≤159KG



DJDX-AL2-1

Aluminum alloy four-fold stretcher


- Product size: 2000*500*150mm
- Folding size: 1070*210*110mm
- Material: Aluminum alloy folding stretcher
- Tube wall thickness: 2.5mm, diameter 32mm
- Weight: 4.6kg
- Load-bearing capacity: ≤159KG



DJAL-MDJCS-6

Scoop Stretcher

- Unfolded size: 2100*440*60mm
- Folded size: 1200*440*90mm
- Load bearing: ≤159KG
- Net weight: 7.56KG
- Gross weight: 8.89KG
- Packing box size: 1705*460*70mm, one piece/carton



Computer Inspection Series



From the emergency room to the inpatient department, and from routine inspections to emergency responses, our cart products, with medical-grade protection standards, ensure zero error in data collection. They support quick expansion of multiple devices and are compatible with mainstream medical systems, allowing real-time synchronization and traceability of every test result. Through the power of technology, we're shortening the diagnostic distance and prolonging the hope of life.



CTAB-PSDJJ-5

Emergency trolley

Specification

- net weight: 136KG
- Material: Aluminium, Steel, ABS
- Concealed sub-work desk, a file box,One oxygen cylinder bracket, one cardiopulmonary resuscitation press plate, and one power socket,Two ABS dirt buckets, one 3-liter sharps box
- Size: 750*480*930mm
- Optional: 850*520*950mm



CTAB-PSDJJ-4



CTAB-PSDJJ-6



CTAB-PSDJJ-13

CTAB-PSDJJ-17

Emergency trolley

Specification

- Material: The main body is made of ABS engineering plastic
- Drawer: Modular design of all-plastic ABS material;
- Column: Plastic steel column
- Fence: Stainless steel three-sided fence
- Car body: Plastic steel gold column
- Central control: Central control locking device



CTAB-JHJJ-4
Size:750*480*930mm



CTAB-JHJJ-2
Size:850*500*950mm



CTAB-SRMJJ-1
Size:625*475*950mm

CTAB-JHSY-2

Anesthesia trolley

Specification

- Material: ABS engineering plastic;
- Drawer: All-plastic ABS material
- Castors: 100mm diameter, 4 silent casters
- Size: 750*480*1300mm
- Optional: 850*520*950mm



CTAB-PSDSY-2
Size:750*480*930mm



CTAB-PSDSY-3
Size:625*480*920mm



CTAB-PSDSX-2
Size:750*480*930mm

CTAB-JHMZ-6

Anesthesia trolley

- Size: 750*480*1300mm
- Optional:850*520*950/1320mm
- Material: ABS engineering plastic;
- Drawer: All-plastic ABS material
- Castors: 100mm diameter, 4 silent casters



BLAB-PSD-1

Medical record trolley

- Style: double drawer,
- Optional: 20/40, 25/50, 30/60 compartments
- Casters: 2 with brakes, 2 without brakes



CTAB-JHMZ-3
Size:750*480*930/1300mm



CTAB-JHMZ-5
Size:810*520*1020mm



CTAB-PSDMZ-5
Size:810*520*1320mm



CTSS-JHWW-1
Size:850*570*880mm



CTSS-MZWW-1
Size:550*550*860mm



CTSS-MZWW-5
Size:900*450*850mm

CTSS-MZ2-9

Stainless steel cart

- Size: 670*450*900mm;
- Castors: 3-inch ABS casters, 2 with brakes
- Frame material: stainless steel 25*0.8mm
- Guardrail material: stainless steel 13*0.8mm
- Panel: 0.8mm or more
- Others: with a stainless steel bowl and a plastic bucket



CTSS-MZ3-3

Stainless steel cart

- Size: 600*400*860mm/
660*440*860mm/
800*480*860mm
- Casters: 4 casters
- Material: Stainless steel
- Guardrail: Stainless steel guardrail

CTSS-MZ2-1

Stainless steel cart

- Size: 600*400*860mm/
660*440*860mm/
800*480*860mm
- Casters: Four 3-inch grey rubber wheels, 2 with brakes
- Material: stainless steel; 0.6mm thickness



CTSS-MZ2-6

Stainless steel cart

- Size: 600*400*860mm/
660*440*860mm/
800*480*860mm
- Casters: four 3-inch grey rubber wheels,
2 with brakes
- Material: stainless steel; 0.6mm thickness

CTSS-MZ2-2

Stainless steel cart

- Size: 600*400*860mm/
660*440*860mm/
800*480*860mm
- Casters: four 3-inch grey rubber wheels, 2 with brakes
- Material: 430 stainless steel; 0.6mm thickness



CTSS-KXD-4

Stainless steel mayo trolley

· Size: 600*400*800-1100mm



CTSS-YKX-5

Stainless steel medical table

· Size: 1000*500*800mm



CTSS-YKX-6

Stainless steel medicine cabinet

· Size: 1500*600*2000mm



XSSS-SDTH-1

Stainless steel sink

· Size: 1800*600*1800mm



CTAB-QB2-1

ABS cart

· Size: 539*366*830mm
640*456*830mm



CTAB-PSD2-2

ABS cart

· Size: 539*366*830mm
640*456*830mm



CTAB-PSD2-3

ABS cart

· Size: 539*366*830mm
640*456*830mm



CTAB-QB3-1

ABS cart

· Size: 539*366*880mm
640*456*900mm



CTAB-PSD3-4

ABS cart

· Size: 539*366*900mm
640*456*900mm



CTAB-PSD3-5

ABS cart

· Size: 539*366*900mm
640*456*900mm



CTAB-PSD3-6

ABS cart

· Size: 539*366*880mm
539*366*900mm



MRSL-KL-1

Elderly recliner

Specification

· Product size: 1100*660*1180mm · Packaging size: 950*700*720mm · Net weight: 41kg

Description: The fabric is PVC calendered imitation leather, steel pipe sprayed, with three adjustable angles, with lockable dining table and 5-inch casters.



MRSL-KL-3

· Product Size: 1400*7300*1180mm



MRSL-KL-2

· Product size: 1100*770*1280mm



MRSL-KL-5

Electric dialysis chair

Specification

- Product size: 900x880x1320mm
- Net weight: 60kg
- Packaging size: 1510x800x610mm
- Trendelenburg Postion(-12°)



MRSL-KL-4

Manual Dialysis Chair

Specification

- Product size: 900x880x1320mm
- Packaging size: 1510x800x610mm
- Net weight: 55kg



MRSL-KL-6

Electric dialysis chair

Specification

- Product size: 1060x870x1300mm
- Cushion height: 560-820mm
- Packaging size: 1610x820x600mm
- Net weight: 60kg



PHG-ZK-1



Accompanying chair

- Material: thickened steel pipe, thickened leather, sponge
- Size: 690*600*1000mm
- Function: can lie down, can sit



PHG-PL-1



Accompanying chair

- Chair state (folded size): 790*630*360mm
- Backrest height: 840mm
- Bed state (opened size): 1880*630*490mm
- Seat thickness: 30mm

