



GRACE MEDICAL INSTRUMENT CO.,LTD

TEL:+86 193 3372 3988

E-mail:supplier@gracemedy.com

Website:www.gracemedy.com

Address:No. 10 Xiwaihuan Road, Zaoqiang County
Hengshui City, Hebei Province, China



Official website



WhatsApp

MEDICAL DEVICE MANUFACTURE
ONE-STOP SERVICE

GRACE MEDY PRODUCT BROCHURES

Medical Accessories

GRACE MEDY

Embrace The Future

Gracemedy Is Here

Escort Your Health and Spread Love



CONTENT

01	GRACE STRENGTH
03	HEADBOARD & FOOTBOARD
06	GUARDRAIL


13	CRANK
14	DINING TABLE
15	CASTERS

17	AIR MATTRESS
19	IV POLE
21	HOSPITAL BED MOTOR

23	HAND CONTROLLER
----	-----------------



GRACE STRENGTH



Grace Medy is a professional manufacturer and designer of hospital beds and accessories, wheelchair and relate medical products. Our factory established in 2009, after years of development, we have created a professional R&D team focusing on electric ICU beds, home care beds, electric wheelchair and all the related accessories. Our comprehensive and intelligent production lines provide you with modern medical and nursing combined products.

Our company has passed CE certification, FDA, and ISO13485 quality management system certification. These certificates are a recognition of our quality.

Our commitment: First-class products, First-class quality, First-class service



HEADBOARD & FOOTBOARD

MATERIAL: ABS/PP/COMPOSITE



CTP-JC-1

- Size:960*520mm
- Material: PP
- Installation distance: 720mm



CTAB-WK-1

- Size:940*580/510mm
- Material: ABS
- Installation distance: 600mm



CTP-CF-7

- Size:960*550/490mm
- Material: PP



CTP-DH-3

- Size:950*540/440mm
- Material: PP
- Installation distance: 900mm



CTAB-WK-05

- Size:957*535/445mm
- Material: ABS
- Installation distance: 900mm



CTFH-CF-3

- Size:940*500mm
- Material:ABS composite stainless steel pipe
- Installation distance: 905-910mm



CTAB-WK-03

- Size:1050*590/530mm
- Material: ABS
- Installation distance: 620mm



CTAB-WK-10

- Size:900*570/520mm
- Material: ABS
- Installation distance: 470mm



CTP-YK-1

- Size:960*550/490mm
- Material: PP



HEADBOARD & FOOTBOARD

CTABWK-08

- Size:940*530mm
- Material: ABS
- Installation distance: 600mm



MORE STYLES



CTP-WK-3



CTP-CF-10



CTFH-CF-1



CTP-CF-2



CTP-CF-9



CTFH-CF-6

GUARDRAIL

MATERIAL: ABS/PP/COMPOSITE



HLAB-HDSX-1

- Size:720 *300mm
- Material: PP



HLP-KKSX-5

- Size:720 *300mm
- Material: PP



HLP-RCZSX-1

- Size:810*356mm
- Material: PP



HLP-JSSX-1

- Size: 710*350*mm
- Material: PP



HLAB-WKLD-1

- Size: 1070*280mm
- Material: PP



HLP-KKLD-1

- Size: 1070*300mm
- Material: PP

GUARDRAIL

MORE STYLES



HLP-CFLD-1



HLP-CFLD-2



HLP-CFSX-2



HLP-CFSX-7



MORE STYLES




HLP-CFSX-8



HLP-CFSX-9



HLP-KKSX-6



HLAL-CFSX-1



HLP-KKSX-3



HLP-CFSX-10

RENDERING






HLAL-HCC4-1

- Size:Top tube length 92cm, down tube length 90cm
- Material: Aluminum alloy



HLAL-HCC5-7

- Size:Top tube length 116.5cm, down tube length 113cm
- Material: Aluminum alloy



HLAL-ZCC6-2

- Size:Top tube length 141.5cm, down tube length 133cm
- Material: Aluminum alloy



HLAL-HCC6-7

- Size:Top tube length 141.5cm, down tube length 133cm
- Material: Aluminum alloy



HLAL-OBC6-1

- Size:Top tube length 141.5cm, down tube length 133cm
- Material: Plastic, aluminum alloy

GUARDRAIL



HLSS-DHC5-2




HLAL-DFC6-17



HLSS-DFHS-1



HLSS-QYHS-1



HLSS-DFLD-1



HLSS-XS15-1

CRANK

MATERIAL: ABS/STAINLESS STEEL



YBAB-FLBT-1

- Material: ABS
- Customizable



YBAB-FLST-1

- Material: ABS
- Customizable



YBSS-LRBT-1

- Material: Stainless steel
- Customizable



YBSS-FLST-1

- Material: Stainless steel
- Customizable

DINING TABLE

MATERIAL: ABS/CARBON STEEL



CZP-LYCB-1

- Size: 880-1060*320*25mm
- Material: ABS
- Features: retractable, waterproof, with hook

CZAB-YKCB-2

- Size: 880-1060*320*25mm
- Material: ABS
- Features: retractable, waterproof, with hook

CZMZ-WKYD-1

- Desktop size: 900*400mm
- Height : 750-1140mm
- Main body material: carbon steel
- Load bearing: 30kg



CZAB-WKYD-1

- Desktop size: 820*450mm
- Height : 840-1080mm
- Main body material: carbon steel
- Load bearing: 30kg

CZMZ-JLYD-1

- Desktop size: 765*380mm
- Height 680-1000mm
- Main body material: carbon steel
- Load bearing: 20kg

CZMZ-JLYD-2

- Desktop size: 600*400mm
- Height 650-1000mm
- Main body material: carbon steel
- Load bearing: 20kg

CASTERS

MATERIAL: POLYAMIDE/STAINLESS STEEL



JL65-YHPE-1

- Application Spoce:Applicable to the treatment vehicles, instruments, hospital beds,etc;
- High strength Polyamide bracket,specious and antispetic; Needle bearing with steel ball loading structure design,
- flexible rotation and good carrying capacity;
- The color of Bracket and wheel can be customized;
- Installation method and dimension can be customized.



JL100-YKQS-1

- Applicable to the treatment vehicles,hospital instruments,baby bed; High strength polyamide frosted bracket, goodappearance,
- antiseptic and rust prevention; Precision ball bearing structure for rotating parts,flexible
- rotation,silent and good loading capacity;
- Directional function can be customized;
- Customizable stainless steel hardware,
- non-magneticcaster can be customized;
- The color of Bracket and wheel can be customized;
- Installation method and dimension can be customized.



JL125-PTTPA-6

- Applicable to hospital beds,nursing beds, emergency vehicles, ect.
- The wheels have built-in precision ball bearings.
- Integrated wheel rack.
- Nylon main drive and wheel core are used, with strong bearing capacity.
- TPU wear-resistant wheel surface, noise reduction and shock absorption.



JL125-PTTPR-11

Customizable



JL75-YHTPR-1

Customizable



JL75-PTTPR-7

Customizable



JL125-PTTPR-15

Customizable





AIR MATTRESS



Anti-bedsores air mattress

CDCS-FLM08-1

- size: can be customized according to needs
- Material: Oxford cloth, PU, composite PU optional
- Thickness: can be customized according to needs
- Lining: brown silk, sponge, sponge + palm silk



CDCS-JY75-4

- Size:1950*900mm
- Thickness:6cm,8cm,10cm,12cm,13cm,15cm
- Fabric: PU waterproof fabric
- Padding: palm fiber+sponge, polyester, high density sponge, high resilience memory foam
- Function: 1 function, 2 functions, toilet hole, optional
- Customization: Thickness, Fabric, Size, Filling

IV POLE

MATERIAL: ABS/STAINLESS STEEL/ALUMINUM ALLOY/CAST IRON



- SLCT-XS-4**

Packing size : 115*35*15cm
Material : Stainless steel
Base : Cast iron,ABS
On the tube : Phi 19 mm
The tube : Phi 25 mm
Adjust way : Knob
Height : 2 meters
Hook : 4
Packing : 8 / box
Weight : 25 kg
Gross weight : 27.3 kg
- SLCT-XS-6**

Packing size : 115*35*15cm
Material : Stainless steel
Base : Cast iron, ABS
On the tube : Phi 19 mm
The tube : Phi 25 mm
Adjust way : Rotation
Height : 2 meters
Hook : 4
Packing : 8 / box
Weight : 22 kg
Gross weight : 24.4 kg
- SLCT-ZJ-5**

Packing size : 123*54*22cm
Material : Stainless steel
Base : Stainless steel pipe
On the tube : Phi 16 mm
The tube : Phi 32 mm
Adjust way : Rotation
Height : 2 meters
Hook : 4
Packing : 5 / box
Weight : 18.3kg
Gross weight : 24kg

- SLCT-XS-7**

Packing size : 93*54*25cm
Material : Stainless steel
Base : Stainless steel pipe,ABS plastic
On the tube : Phi 19 mm
In the tube : Phi 25 mm
The tube : Phi 32 mm
Adjust way : Rotation
Height : 2 meters
Hook : 4
Packing : 5/ box
Weight : 21.8 kg
Gross weight : 26.2kg
- SLCT-XS-8**

Packing size : 122*40*22cm
Material : Stainless steel
Base : Cast iron, mat
On the tube : Phi 19 mm
Phi 25 mm
The tube : Phi 32 mm
Phi 40 mm
Adjust way : Rotation
Tray : Aluminum alloy
Accessories :
The storage basket : 1
Mute the brake wheel : 5
Height : 1.5 M to 2.4 M
Hook: 4

- SLCT-XS-20**

Tabletop : Aluminum alloy
Size : 255*145mm
Height from ground : 90cm
Pole : stainless steel
Basket : Iron wire sprayed plastic
Size : 28*16*11cm
Casters : 2.5-inch silent casters
Base : ABS plastic, inner cast iron plate

IV POLE

HOSPITAL BED MOTOR



- Protection level: Max. IP66
- Noise level: ≤ 50dB
- Other options: Mechanical anti-pinch / Hall sensor
- Optional hand crank / manual release

Speed and thrust comparison table (2600rpm)

Load (N)			Extension current (A)		Extension speed (mm/s)		Lead	Number of threads
Thrust	Pull	Self-locking force	No load	Load	No load	Load		
8000	4000	6000	0.4	4.2	4.3	3.6	3	1
6000	4000	6000	0.3	3.9	5.8	4	4	1
4000	2000	2400	0.4	1.9	8.2	7.3	6	1
2000	1000	800	0.4	1.5	15.7	14.3	6	2
1000	1000	400	0.4	1.2	30.2	27.3	12	2

DJCT-JC02-2

- Max.load in push &self-locking(N):6000
- Max.static bending moment(Nm):1000
- Max.dynamic bending moment(Nm):500
- Max.speed(mm/s):20
- Installation dimensions:
L=S+200(S<400)
L=S+220(S≥400)
- IP grade:IPX4
- Option:
With hall sensor(Apply for synchronize or memory function)



DJCT-JC02-3

- Max.load in push &self-locking(N):2000
- Max.static bending moment(Nm):900
- Max.dynamic bending moment(Nm):500
- Max.speed(mm/s):10
- Installation dimensions:
L=S/2+150(S≤400)
L=S/2+200(S>400)
- IP grade:IPX4
- Option:
With hall sensor(Apply for synchronize or memory function)



DJCT-JS02-1

- Maximum thrust : 6000N*2
- Back stroke: Up to 87mm
- Leg stroke: Up to 87mm
- Olnput voltage: 110/230VAC
- Ingress protection: IPX4
- Mounting distance: 581+3mm

HAND CONTROLLER



