



GRACE MEDICAL INSTRUMENT CO.,LTD

TEL:+86 193 3372 3988

E-mail:supplier@gracemedy.com

Website:www.gracemedy.com

Address:No. 10 Xiwaihuan Road, Zaoqiang County
Hengshui City, Hebei Province, China



Official website



WhatsApp

MEDICAL DEVICE MANUFACTURE
ONE-STOP SERVICE

GRACE MEDY PRODUCT BROCHURES

Rehabilitation Equipment

GRACE MEDY

Embrace The Future

Gracemedy Is Here

Escort Your Health and Spread Love



CONTENT

01	GRACE STRENGTH
03	ELECTRIC WHEELCHAIR
15	MANUAL WHEELCHAIR

29	COMMODE CHAIR
32	SHOWER CHAIR
33	WALKING AID

35	SHOPPING CART
36	WALKING STICK
40	ARMPIT CRUTCH



GRACE STRENGTH



Grace Medy is a professional manufacturer and designer of hospital beds and accessories, wheelchair and relate medical products. Our factory established in 2009, after years of development, we have created a professional R&D team focusing on electric ICU beds, home care beds, electric wheelchair and all the related accessories. Our comprehensive and intelligent production lines provide you with modern medical and nursing combined products.

Our company has passed CE certification, FDA, and ISO13485 quality management system certification. These certificates are a recognition of our quality.

Our commitment: First-class products, First-class quality, First-class service



LYPS-YWRECG-4

Electric wheelchair



SPECIFICATION	
SIZE	945*640*875mm
Weight	55kg
Load-Bearing	130kg

Front/rear wheels: 6-inch swivel wheels/10-inch solid wheels

Braking system: EPS automatic braking system

Obstacle height: 30mm

Seat lift stroke: 180mm

Handrail flip angle: 0-90°

Battery/range: Lithium battery 10Ah/8km

Charging time: 5h

500W Wanshuo motor Motor: 350w*2

Gears: Slow 5th / Fast 5th

Maximum speed: 6.3km/h

Boot mode: mechanical key unlocking

Lights: Headlights

LYPS-ZYECG-3

Electric wheelchair



SPECIFICATION	
SIZE	1120*700*1060mm
Weight	113kg
Load-Bearing	100kg

Seat width 44cm; The backrest is 48cm high

Front/rear: 10 inches/16 inches

Lead-acid battery: 24v 58Ah

Charging time: 5-8h

Endurance: 35km

Turbo pot rod motor: 500W*2

Speed: 1-12km/h (5 gears adjustable)

Material: carbon steel skeleton

Wooden box packing: 1 car/box

Packing size: 98*75*99cm

SPECIFICATIONS	Overall Length cm	Overall Height cm	Overall Width cm	Back Height cm	Seat Height cm	Seat Depth cm	Seat Width cm	Folded Width cm	Loading kgs	G.W. kgs	Packing Size cm
LYPS-YWRECG-4	94.5	87.5	64	32.5	55	41	47	*	130	55	99*68*58
LYPS-ZYECG-3	112	106	70	48	58	45	44	*	100	113	98*75*99

LYPS-ZYECG-5

Electric wheelchair



SPECIFICATION	
SIZE	1120*700*1140mm
Weight	118kg
Load-Bearing	100kg

Seat width 48cm; The backrest is 63 cm high

Seat height from the ground 59cm; Seat depth 47cm

Material: carbon steel frame

Front/rear: 10-inch pneumatic tyres/16-inch pneumatic tyres

Turbo worm motor: 500W*2,
electromagnetic brake

Speed: 1-12km/h (5 gears adjustable)

Lead-acid battery: 24v 58Ah

Endurance: 35-50km

Charging time: 5-8h

Packed in wooden box, 1 set/box

Packing size: 89*72*81cm (excluding tray)

LYPS-ZYECG-5

LYPS-ZYECG-4

Electric wheelchair



SPECIFICATION	
SIZE	1260*680*1190mm
Weight	126kg
Load-Bearing	100kg

Seat width 50cm; The backrest is 82 cm high

Seat height from the ground 58cm; Seat depth 50cm

Material: carbon steel frame

Front/rear: 10-inch pneumatic tyres/16-inch pneumatic tyres

Turbo worm motor: 500W*2

Speed: 1-12km/h (5 gears adjustable)

Lithium battery: 24v 58Ah

Endurance: 35km

Charging time: 5-8h

Packed in wooden box, 1 set/box

Packing size: 98*75*99cm

SPECIFICATIONS	Overall Length cm	Overall Height cm	Overall Width cm	Back Height cm	Seat Height cm	Seat Depth cm	Seat Width cm	Folded Width cm	Loading kgs	G.W. kgs	Packing Size cm
LYPS-ZYECG-5	112	114	70	63	58	45	48	*	100	118	98*75*99
LYPS-ZYECG-4	126	119	68	65	58	45	48	*	100	126	98*75*99

LYPS-YFECG-1

Electric wheelchair



LYPS-YFECG-1

SPECIFICATION

SIZE	975*630*1080mm
Weight	129kg
Load-Bearing	100kg

Front/rear wheels: 10-inch omni-directional wheels/12-inch solid wheels

Seat lift stroke: 180mm

Obstacle height: 80mm

Handrail flip angle: 0-90°

Motor: 450W*2

Gears: Slow 5 gears / Fast 5 gears

Maximum speed: 7km/h

Lithium battery: 70Ah/36km

Braking system: EPS automatic braking system

Boot mode: password/fingerprint

Lights: 4-level high-brightness headlights/steering taillights

LYPS-SWDB-2

Electric wheelchair



SPECIFICATION

SIZE	1280*650*1300mm
Weight	120kg
Load-Bearing	130kg

Battery: 40AH lead-acid / 24V

Motor: 350W

Packaging size (length width height): 1340*720*1180

Speed: 12km/h

Endurance: 20-25km

Front and rear 10-inch inflatable tires

SPECIFICATIONS	Overall Length cm	Overall Height cm	Overall Width cm	Back Height cm	Seat Height cm	Seat Depth cm	Seat Width cm	Folded Width cm	Loading kgs	G.W. kgs	Packing Size cm
LYPS-YFECG-1	97.5	108	63	30~50	50.5~68.5	45	40	*	100	129	95*70*95
LYPS-SWDB-2	128	130	65	65	70	45	46	*	130	120	134*72*118

LYAL-ZYECG-2

Electric wheelchair

SIZE	MATERIAL
1020*600*980mm	Aviation aluminum alloy
SPECIFICATIONS	
7 inch magnesium alloy solid wheel; rear wheels 11-inch magnesium alloy solid wheels; Brushless motor 150W*2;Lithi-um battery 24V 6AH; 6 km/h; Endurance 10-30 km; Net weight: 16kg; Capacity: 100kg	



LYAL-ZYECG-6

Electric wheelchair

SIZE	MATERIAL
1100*630*960mm	Aviation aluminum alloy
SPECIFICATIONS	
8-inch solid front wheel; 12-inch inflatable rear wheels; 12AH lithium battery; Domestic brush motor: 250W; Speed: 6km/h; Endurance: 15-20km; electromagnetic brakes; Maximum load: 120kg; Weight: 27.8kg	



LYPS-SWCECG-7

Electric wheelchair

SIZE	MATERIAL
960*600*980mm	Aviation aluminum alloy
SPECIFICATIONS	
7-inch front wheels PU foam + 9-inch rear wheels; PU foam foot pedal; Dual-motor drive (brushless mechanical motor); Motor: 24V/180W; Maximum speed: 6km/h	



LYAL-YWRECG-1

Electric wheelchair

SIZE	MATERIAL
1020*600*935mm	Aviation aluminum alloy
SPECIFICATIONS	
8-inch solid front wheel; 12-inch inflatable rear wheels; Speed: ≤6km/h; Driving distance: ≥15km; electronic standing slope;12AH lithium battery; electromagnetic brakes; Motor: 24V 250W; Net weight: 24.9kg	

SPECIFICATIONS	Overall Length cm	Overall Height cm	Overall Width cm	Back Height cm	Seat Height cm	Seat Depth cm	Seat Width cm	Folded Width cm	Loading kgs	G.W. kgs	Packing Size cm
LYAL-ZYECG-2	102	98	60	44	50	42	45	32	100	20	61*36*78.5
LYPS-SWCECG-7	96	98	60	52	53	42	41	36	100	18	56*34*81

SPECIFICATIONS	Overall Length cm	Overall Height cm	Overall Width cm	Back Height cm	Seat Height cm	Seat Depth cm	Seat Width cm	Folded Width cm	Loading kgs	G.W. kgs	Packing Size cm
LYAL-ZYECG-6	110	96	63	32	55	45	49	38	100	31	65*42*79
LYAL-YWRECG-1	102	93.5	60	45	48	45	44	38	100	28	65*42*79

LYAL-SWCECG-5

Electric wheelchair

SIZE	MATERIAL
1100*630*960mm	Aviation aluminum alloy
SPECIFICATIONS	
8-inch front wheels; 12-inch rear wheels; 24V12AH lithium battery; Charging time: 5-6H; 250W*2 brushed motor; Speed: 1-6km per hour (5 gears adjustable); Endurance: 15km; Load: 100kg; Net weight: 24.5kg+3kg (lithium battery)	



LYAL-SWCECG-7

Electric wheelchair

SIZE	MATERIAL
1100*650*1250mm	Aviation aluminum alloy
SPECIFICATIONS	
8-inch aluminum alloy solid wheel front wheel size; 12-inch aluminum alloy pneumatic wheel rear wheel size; 12AV20A lithium battery (24V); Charging time: 6-8h; electric recline; 250W*2 (brushed) electromagnetic brake motor; Cruising range: 15-25km; Load: 100kg; Net weight: 32kg	



LYPS-JYECG-11

Electric wheelchair

SIZE	MATERIAL
1130*630*1280·mm	Carbon steel
SPECIFICATIONS	
10-inch front wheels; 16-inch rear wheels; 24V12AH lead-acid battery; Charging time: 5-6h; 250W*2 brushed motor; SPEED: 1-6KM (5 gears adjustable); Cruising range: 15km; electric leg lift; Capacity: 100kg	



LYPS-YFECG-3

Electric wheelchair

SIZE	MATERIAL
995*650*890mm	Aviation aluminum alloy
SPECIFICATIONS	
8-inch swivel wheel front wheel; 20Ah lead-acid battery; Endurance: 30km; 350W*20,000 Shuo motor; Maximum speed: 7km/h; EPS automatic braking system; Charging time: 5h; Vehicle weight: 48kg; The vehicle bears 100kg	

SPECIFICATIONS	Overall Length cm	Overall Height cm	Overall Width cm	Back Height cm	Seat Height cm	Seat Depth cm	Seat Width cm	Folded Width cm	Loading kgs	G.W. kgs	Packing Size cm
LYAL-SWCECG-5	110	96	63	32	55	45	49	38	100	31	65*42*79
LYPS-JYECG-11	113	128	63	78	55	46	44	37	100	50	81*40*83

SPECIFICATIONS	Overall Length cm	Overall Height cm	Overall Width cm	Back Height cm	Seat Height cm	Seat Depth cm	Seat Width cm	Folded Width cm	Loading kgs	G.W. kgs	Packing Size cm
LYAL-SWCECG-7	110	125	65	40	53	45	49	42	100	36	65*46*80
LYPS-YFECG-3	99.5	89	65	50	52.5	45	40	*	100	55	99*68*68



LYPS-ZYECG-8

Electric wheelchair

- Material: Carbon steel metal paint
- Front wheels: 7-inch solid tire wheels
- Rear wheels: 12-inch solid tire wheels
- Battery:Lead-acid battery 24V /12AH
- Motor: 250W*2pcs
- Endurance: 15-20KM
- Speed: 1-6km/h (5 gears adjustable)
- Charging time: 5-6h



LYI-SWCECG-1

Electric wheelchair

- Material: Carbon steel metal paint
- Front wheels: 10-inch solid rubber tire wheels
- Rear wheels: 24-inch solid tire wheels
- Battery: Lead-acid 24V/12AH
- Controller: 360° intelligent gimbal controller
- Motor: 250W*2pcs
- Endurance: 15-25KM
- Speed: 1-6km/h (5 gears adjustable)
- Brakes: Electronic brakes
- Charging time: 6-8h



LYPS-SWCECG-3

Electric wheelchair

- Material: Carbon steel metal paint
- Front wheels: 10-inch solid rubber tire wheels
- Rear wheels: 16-inch inflatable tire wheels
- Battery: Lead-acid 24V/12AH
- Controller: 360° intelligent gimbal controller
- Motor: 250W*2pcs
- Endurance: 15-25KM
- Speed: 1-6km/h (5 gears adjustable)
- Brakes: Electronic brakes
- Charging time: 6-8h

SPECIFICATIONS	Overall Length cm	Overall Height cm	Overall Width cm	Back Height cm	Seat Height cm	Seat Depth cm	Seat Width cm	Folded Width cm	Loading kgs	G.W. kgs	Packing Size cm
LYPS-ZYECG-8	106	94	60	44	50	42	45	36	100	32	74*40*43
LYI-SWCECG-1	125	98	68	45	53	46	47	38	100	47	80*39*82
LYPS-SWCECG-3	111	96	65	45	53	46	45	43	100	47	80*39*81



LYPS-JYECG-8

Electric wheelchair

- Controller: 360° intelligent gimbal controller
- Motor: 250W*2pcs
- Battery: 24V12AH(20AH)
- Voltage: 220-240V
- Endurance: 15-20KM
- Top speed: 6km/h
- Brakes: Fully automatic electronic brakes
- Rear wheel size: 16-inch pneumatic tires, aluminum alloy wheels
- Front wheel size: 10-inch solid tires, aluminum alloy wheels



LYPS-FLCECG-5

Electric wheelchair

- Maximum speed: 1-6km/h
- Maximum driving incline: 13°
- Endurance: 15-25km
- Braking system: Intelligent electronic brake
- Battery: Lead-acid battery 24V 12AH/20AH*2pcs
- Motor: 250W*2pcs
- Rear wheel size: 16-inch pneumatic tires wheels
- Front wheel size: 10-inch solid tires wheels



LYPS-YWRECG-1


Electric wheelchair

- Material: carbon steel frame, sprayed surface treatment
- Speed: 1-6km/h; Driving distance: ≥15km
- Front wheels: 10-inch white solid wheels
- Rear wheels: 22-inch white pneumatic tyres
- Hand electric clutch switching
- Brushed motor: 24V 250w
- Lead-acid battery: 12AH
- Red or blue seat cushion back cushion

SPECIFICATIONS	Overall Length cm	Overall Height cm	Overall Width cm	Back Height cm	Seat Height cm	Seat Depth cm	Seat Width cm	Folded Width cm	Loading kgs	G.W. kgs	Packing Size cm
LYPS-JYECG-8	128	98	74	45	53	46	56	38	100	47.38	86*43*82
LYPS-FLCECG-5	111	92	65	45	47	43	43	37	100	40.3	79*40*79
LYPS-YWRECG-1	117	96	66	48	48	45.5	44.5	44	100	46.7	84.5*46.5*79

LYAL-TKMCG-8

Manual wheelchair



SPECIFICATION	
SIZE	1070*540*1000mm
Weight	24.6kg
Load-Bearing	100kg

Material: Aluminum alloy frame

Seat back cushion: split sponge seat, adjustable headrest with shoulder straps

Backrest: The wheelchair can adjust the reclining angle (90 ° -160 °)

and the sitting angle (90 ° -45 °). The backrest is equipped with adjustable fixed pads on both sides

LYPS-M KNT-1

LYAL-TKMCG-7

Manual wheelchair



SPECIFICATION	
SIZE	1035*620*940mm
Weight	15.8kg
Load-Bearing	100kg

Material: aluminum alloy frame, surface spraying

Seat back cushion: 5cm thickened high resilience decompression seat cushion

Handrails: Fixed handrails

Footrest: Fixed footrest, thickened plastic pedal

Safety configuration: with seat belts

Brake: toggle brake

The front wheels are made of 8-inch PVC material, and the rear wheels are 24-inch inflatable rear wheels, and the plastic hand ring is raised

SPECIFICATIONS	Overall Length cm	Overall Height cm	Overall Width cm	Back Height cm	Seat Height cm	Seat Depth cm	Seat Width cm	Folded Width cm	Loading kgs	G.W. kgs	Packing Size cm
LYAL-TKMCG-8	103	92	49	46	46	46	37	29	100	24.4	94*29*85
LYAL-TKMCG-7	103.5	94	62	41	41	55	43	28	100	15.8	103*27*98

LYAL-TKMCG-2

Manual wheelchair

SIZE	LOAD-BEARING
1095*960*635mm	100kg

SPECIFICATIONS
Aluminum alloy frame
8-inch PU wheels in the front and 24-inch pneumatic tires in the rear
The push support is detachable, and the handrail can be lifted and lifted
PU seat cushion with thickened removable sponge cushion



LYAL-TKMCG-4

Manual wheelchair

SIZE	LOAD-BEARING
960*470*895mm	75kg

SPECIFICATIONS
Armrest: The armrest can be turned back, and the height can be adjusted in 5 gears
Seat belts: Pelvic seat belts are adjustable
Leg protection straps: Comes with shin straps
The front wheels are 6-inch solid tires



LYAL-TKMCG-9

Manual wheelchair

SIZE	LOAD-BEARING
1070*580*920mm	100kg

SPECIFICATIONS
Handrails: Flip-back handrails
Seat back cushion: double mesh seat back cushion is detachable, good air permeability
Backrest: The backrest is foldable
Footrest: detachable movable pedal, plastic pedal



LYAL-MKNT-1

Manual wheelchair

SIZE	LOAD-BEARING
1035*620*940mm	100kg

SPECIFICATIONS
Frame material: Aluminum alloy material
Frame color: blue, yellow
Adjustable height armrest
Pedal: Pedal, foldable

SPECIFICATIONS	Overall Length cm	Overall Height cm	Overall Width cm	Back Height cm	Seat Height cm	Seat Depth cm	Seat Width cm	Folded Width cm	Loading kgs	G.W. kgs	Packing Size cm
LYAL-TKMCG-2	109.5	96	63.5	40	52	43	46	28	100	20	83*28*98
LYAL-TKMCG-9	107	92	58	41	50	43	41	26	100	17.9	83*27*77

SPECIFICATIONS	Overall Length cm	Overall Height cm	Overall Width cm	Back Height cm	Seat Height cm	Seat Depth cm	Seat Width cm	Folded Width cm	Loading kgs	G.W. kgs	Packing Size cm
LYAL-TKMCG-4	96	89.5	47	36	53	38	35	28.5	75	15.5	75*29*92
LYAL-MKNT-1	103.5	94	62	46	46	37	37	29	75	22.8	91*30*96

LYAL-TKMCG-3

 Manual wheelchair

SIZE	LOAD-BEARING
1030*675*940mm	100kg
SPECIFICATIONS	
Seat Belt: Pelvic Seat Belt	
Front wheel: aluminum alloy fork, 8-inch PU solid tire	
Rear wheel: 24 inch inflatable tire for rear wheel	



LYAL-TKMCG-5

 Manual wheelchair

SIZE	LOAD-BEARING
750*480*940mm	100kg
SPECIFICATIONS	
Material: Aluminum alloy frame	
Seat back cushion: double-layer seat back cushion	
Front and rear wheels: 4-inch solid universal wheels for the front wheel, 20 inch gray inflatable tires for the rear wheel, with anti rollover and obstacle crossing functions	



LYAL-TKMCG-20

 Manual wheelchair

SIZE	LOAD-BEARING
1050*630*930mm	100kg
SPECIFICATIONS	
Material: Aluminum alloy frame	
Seat back cushion: black flame-retardant Oxford cloth, detachable and thickened 5cm back cushion	
Rear wheel: 24 inch inflatable tire for rear wheel	



LYPS-MKNT-1

 Manual wheelchair

SIZE	LOAD-BEARING
940*290*850mm	100kg
SPECIFICATIONS	
Material: Thickened steel	
Seat cushion: It is a hydraulically adjustable angle bee net cushion with height-adjustable headrest; height-adjustable handrails;	
Guards: Beenet side baffles on both sides	

SPECIFICATIONS	Overall Length cm	Overall Height cm	Overall Width cm	Back Height cm	Seat Height cm	Seat Depth cm	Seat Width cm	Folded Width cm	Loading kgs	G.W. kgs	Packing Size cm
LYAL-TKMCG-3	103	94	67.5	43	55	45	50	26	100	17.7	103*27*82
LYAL-TKMCG-20	105	93	63	42	54	44	46	27	100	18	83*28*78

SPECIFICATIONS	Overall Length cm	Overall Height cm	Overall Width cm	Back Height cm	Seat Height cm	Seat Depth cm	Seat Width cm	Folded Width cm	Loading kgs	G.W. kgs	Packing Size cm
LYAL-TKMCG-5	75	94	48	44	45	34	30	34	100	15.5	72*35*64
LYPS-MKNT-1	103	92	49	46	46	46	37	29	100	24.4	94*29*85

LYT-TKMCG-1

Manual wheelchair

SIZE	LOAD-BEARING
1105*665*1240mm	100kg
SPECIFICATIONS	
Material: Carbon Fiber	
Carbon fiber footrest with adjustable height	
4-inch 5-inch solid core small front wheel optional, inflatable rear wheel 24/26 inch optional	



LYAL-KHCG-1

Manual wheelchair

SIZE	LOAD-BEARING
1020*950*550mm	115kg
SPECIFICATIONS	
Material: Aluminum alloy	
Wheels: Both front and rear wheels are 8-inch wheels	
Net weight: 9kg	
Color: Grey	



LYPS-MKCG-4

Manual wheelchair

SIZE	LOAD-BEARING
800*520*900mm	100kg
SPECIFICATIONS	
Aluminum alloy frame	
Folding size 80 * 28 * 52	
Front 6-inch and rear 12 inch honeycomb elastic wheels	
Breathable mesh cushion	



LYPS-XRCG-1

Manual wheelchair

SIZE	LOAD-BEARING
1100*920*750mm	150kg
SPECIFICATIONS	
Material: welded steel pipe	
Front wheel: 8-inch solid front wheel	
Rear wheel: 24 inch solid tire nylon rear wheel	
Armrest: Plastic armrest	

SPECIFICATIONS	Overall Length cm	Overall Height cm	Overall Width cm	Back Height cm	Seat Height cm	Seat Depth cm	Seat Width cm	Folded Width cm	Loading kgs	G.W. kgs	Packing Size cm
LYT-TKMCG-1	110.5	124	66.5	39	54	45	44	25	100	24.6	97*32*100
LYPS-MKCG-4	80	90	52	53	45	41	40	28	100	11	56*26*80

SPECIFICATIONS	Overall Length cm	Overall Height cm	Overall Width cm	Back Height cm	Seat Height cm	Seat Depth cm	Seat Width cm	Folded Width cm	Loading kgs	G.W. kgs	Packing Size cm
LYAL-KHCG-1	102	55	95	45	49	42	45.7	25	115	9.8	35*29*82
LYPS-XRCG-1	110	75	92	30	52	42	56	25	150	21.7	92*26*96

LYAL-KHCG-2

Manual wheelchair

SIZE	LOAD-BEARING
1000*950*550mm	115kg
SPECIFICATIONS	
Material: Aluminum alloy	
Wheels: Both front and rear wheels are 8-inch wheels	
Color: Grey	



LYPS-LHMCG-1

Manual wheelchair

SIZE	LOAD-BEARING
800*900*590mm	115kg
SPECIFICATIONS	
6-inch solid front wheels	16 inch solid rear wheels
Mesh seats and backrests	
With folded back	
Hand brake and nurse brake	



LYPS-XRCG-4

Manual wheelchair

SIZE	LOAD-BEARING
1050*620*840mm	100kg
SPECIFICATIONS	
Frame: Silver carbon steel frame	
Front wheel: 8-inch solid rubber tire	
Rear wheel: 24 inch rear wheel	
Brake: Galvanized handbrake	



LYPS-HGMCG-19

Manual wheelchair

SIZE	LOAD-BEARING
880*580*880mm	100kg
SPECIFICATIONS	
Front wheel: 6-inch rubber wheels	
Rear wheel:16-inch engineering plastic integral wheel hub, polyurethane non inflatable tire	
The brake adopts steel parking brake and manual control brake	

SPECIFICATIONS	Overall Length cm	Overall Height cm	Overall Width cm	Back Height cm	Seat Height cm	Seat Depth cm	Seat Width cm	Folded Width cm	Loading kgs	G.W. kgs	Packing Size cm
LYAL-KHCG-2	100	55	95	45	49	42	45.7	28	115	10	33*29*82
LYPS-XRCG-4	105	84	62	43	41	43	45	27	100	15.2	93.5*27.5*89

SPECIFICATIONS	Overall Length cm	Overall Height cm	Overall Width cm	Back Height cm	Seat Height cm	Seat Depth cm	Seat Width cm	Folded Width cm	Loading kgs	G.W. kgs	Packing Size cm
LYPS-LHMCG-1	80	59	90	40	48	41	49	26	100	12.25	72*26*72
LYPS-HGMCG-19	88	88	58	43	45	40	43	26	100	14.3	78*30*78.5

LYPS-NKMCG-10

 Manual wheelchair

SIZE	LOAD-BEARING
1070*890*660mm	100kg
SPECIFICATIONS	
Front wheel: 7-inch rubber wheels	
Rear wheel: 24 inch integrated wheel belt with handrail	
Frame: Black carbon steel spray painted frame	



LYPS-LHMCG-6

 Manual wheelchair

SIZE	LOAD-BEARING
1030*670*880mm	100kg
SPECIFICATIONS	
Front wheel: 7-inch rubber wheels	
Rear wheel: 24 inch spoke wheel gray tire	



LYPS-HGCG-3

 Manual wheelchair

SIZE	LOAD-BEARING
1000*650*880mm	100kg
SPECIFICATIONS	
Front wheel: 7-inch PVC black wheel	
Rear wheel: 24 inch PVC solid tire	
Seat cushion: black Oxford cloth	



LYPS-HGMCG-5

 Manual wheelchair

SIZE	LOAD-BEARING
1030*660*910mm	100kg
SPECIFICATIONS	
Front wheel: 7-inch rubber wheels	
Rear wheel: 24 inch spoke rubber tire with handrail	
Iron pedal	
Material: Steel frame	

SPECIFICATIONS	Overall Length cm	Overall Height cm	Overall Width cm	Back Height cm	Seat Height cm	Seat Depth cm	Seat Width cm	Folded Width cm	Loading kgs	G.W. kgs	Packing Size cm
LYPS-NKMCG-10	107	89	66	43	45	42	45	27	100	15.8	93*27.5*96
LYPS-HGCG-3	100	88	65	43	46	43	45	*	100	15.1	93*28*93

SPECIFICATIONS	Overall Length cm	Overall Height cm	Overall Width cm	Back Height cm	Seat Height cm	Seat Depth cm	Seat Width cm	Folded Width cm	Loading kgs	G.W. kgs	Packing Size cm
LYPS-LHMCG-6	103	88	67	42	42	42	46	27	100	14.95	96*89.5*26
LYPS-HGMCG-5	103	91	66	41	50	41	45	26	100	15.15	93*27.5*96



LYPS-NKMCG-7

Manual wheelchair

- Size:1300*660*880mm
- Load-bearing: 100kg



LYPS-HGMCG-14

Manual wheelchair

- Size:1030*660*880mm
- Load-bearing: 100kg



LYPS-HGMCG-13

Manual wheelchair

- Size:920*630*910mm
- Load-bearing: 100kg



LYDD-NKMCG-4

Manual wheelchair

- Size:1030*660*880mm
- Load-bearing: 100kg



LYPS-NKMCG-11

Manual wheelchair

- Size:1070*890*660mm
- Load-bearing: 100kg



LYDD-NKMCG-1

Manual wheelchair

- Size:1000*680*860mm
- Load-bearing: 100kg



LYDD-HGMCG-8

Manual wheelchair

- Size:1050*620*840mm
- Load-bearing: 100kg



LYPS-HGQT-2

Manual wheelchair

- Size:1020-1620*690*880-1200mm
- Load-bearing: 100kg

SPECIFICATIONS	Overall Length cm	Overall Height cm	Overall Width cm	Back Height cm	Seat Height cm	Seat Depth cm	Seat Width cm	Folded Width cm	Loading kgs	G.W. kgs	Packing Size cm
LYPS-NKMCG-7	103	88	66	40	48	40	45	26	100	16.5	93*28*90
LYPS-HGMCG-14	103	88	66	40	48	40	45	26	100	18	93*27*93
LYPS-HGMCG-13	920	91	63	40	40	40	44	25	100	14.1	78*28*80
LYDD-NKMCG-4	103	88	66	41	50	41	45	25	100	15.15	93.5*27.5*89
LYPS-NKMCG-11	107	89	66	42	47	42	45	26	100	18	93*27*93
LYPS-NKMCG-1	100	86	68	42	42	42	45	26	100	18	93*27*93
LYDD-HGMCG-8	105	84	62	41	43	41	45	26	100	16.3	93.5*27.5*89
LYPS-HGQT-2	162	88	69	74	46	42	45	*	100	21.5	93*27*93



ZBSS-YK-7

- Commode chair
- Size: 45*54*78-88cm
 - Material: Stainless steel



ZBPS-CF-1

- Commode chair
- Size: 85*48*95cm
 - Material: Chrome steel



ZBDD-CF-1

- Commode chair
- Size: 85*48*95cm
 - Material: Carbon steel plating



ZBSS-JYJ-4

- Commode chair
- Size: 40*54*74-84cm
 - Material: Stainless steel



ZBDD-YZ-1

- Commode chair
- Size: 83*50*100-107cm
 - Material: Carbon steel plating



ZBSS-YK-6

- Commode chair
- Size: 53*43*64-74cm
 - Material: Stainless steel



ZBSS-YK-2

- Commode chair
- Size: 53*43*64-74cm
 - Material: Stainless steel



ZBSS-YK-1

- Commode chair
- Size: 52*44*47.5-55cm
 - Material: Stainless steel



ZBDD-JYJ-2

- Commode chair
- Size: 54*45*72-82cm
 - Material: Chrome steel



ZBPS-YK-1

- Commode chair
- Size: 52*46*79cm
 - Material: Carbon steel spraying

POTTY
CHAIRS

ZBSS-ZC-3

- Commode chair
- Height: 38 cm
 - Material: Stainless steel



ZBSS-YK-3

- Commode chair
- Size: 45*54*78-88cm
 - Material: Stainless steel



ZBSS-ZC-4

- Commode chair
- Height: 38 cm
 - Material: Stainless steel



ZBSS-ZC-5

- Commode chair
- Size: 35*32*35cm
 - Material: Stainless steel



ZBAL-BK-1

- Commode chair
- Size: 77*55.5*102-107cm
 - Material: Stainless steel



MYAL-YK-4

- Shower chair



MYAL-YK-3

- Shower chair



MYAL-YK-1

- Shower chair



MYAL-JLJ-1

- Shower chair

POTTY
CHAIRS

SHOWER
CHAIR



ZXSS-YBCR-1

Walking aid

- Height: 72-82cm
- Material: Stainless steel



ZXSS-YBCR-2

Walking aid

- Height: 75-85cm
- Material: Stainless steel



ZXSS-YBCR-3

Walking aid

- Height: 75-85cm
- Material: Stainless steel



ZXAL-YBCR-6

Walking aid

- Height: 73-92cm
- Material: Aluminum alloy



ZXAL-ZCCR-1

Walking aid

- Material: Aluminum alloy



ZXAL-ZCCR-2

Walking aid

- Height: 74-92cm
- Material: Aluminum alloy



ZXAL-TSTRT-1

Walking aid

- Material: Aluminum alloy



ZXAL-YBCR-1

Walking aid

- Height: 73-83cm
- Material: Aluminum alloy



ZXAL-YBCR-2

Walking aid

- Height: 80cm
- Material: Aluminum alloy



ZXAL-YBCR-3

Walking aid

- Height: 75-87cm
- Material: Aluminum alloy



ZXAL-YBCR-5

Walking aid

- Height: 73-92cm
- Material: Aluminum alloy



GWTG-DY-1
● Shopping cart
• Size:52*35*85cm



WGAL-DY-4
● Shopping cart
• Size:60*52*86-92cm



GWTG-MZ-1
● Shopping cart
• Size:60*75*82-92cm



GWTG-MZ-2
● Shopping cart
• Size:60*75*82-92cm



GZSS-MZSZ-4
● Walking stick
• Height: 71-91cm
• Material: Stainless steel



GZSS-MZSZ-2
● 4 Legged walking stick
• Height: 66-88 cm
• Material: Stainless steel



GZSS-MZSZ-3
● 3 Legged walking stick
• Height: 68-88 cm
• Material: Stainless steel



GZAL-MZSZ-13
● Walking stick
• Height: 82-92 cm
• Material: Aluminum alloy



GZAL-MZSZ-11
● Walking stick
• Height: 82-92 cm
• Material: Aluminum alloy



GZAL-MZMZ-1
● 4 Section blind cane
• Height: Folding length 33.5CM, unfolded length 125CM
• Material: Aluminum alloy



GZAL-MZSZ-10
● 4 Legged walking stick



GZAL-MZSZ-15
● Walking stick



GZAL-MZSZ-16
● Walking stick



GZAL-MZSZ-8
● 4 Legged walking stick



CZAL-MZSZ-11
● 4 Legged walking stick



GZSS-MZSZ-1
● 4 Legged walking stick



GZAL-MZSZ-5
● Walking stick
• Height: 75-98cm cm
• Material: Aluminum alloy



GZAL-MZSZ-2



GZAL-MZSZ-7



GZAL-MZSZ-1



GZAL-MZSZ-9

WALKING
AID



GZAL-MZSZ-4

 Walking stick with light

- Height: 80-100 cm
- Material: Aluminum alloy



GZSS-MZSD-1

 Walking stick

- Height: 85 cm
- Material: Stainless steel



GZSS-YBYX-1

 Underarm crutch

- Large size height: 125 cm
- Medium height: 115 cm
- Small height: 105 cm



CZAL-JHYX-1

 Underarm crutch



CZAL-JHYX-2

 Underarm crutch



GZAL-JHYX-4

 Underarm crutch



GZAL-MZZG-3

 Elbow crutch




GZAL-MAZG-4

 Elbow crutch



GZAL-DRCR-3

 Elbow crutch



GZAL-DSZG-1

 Elbow crutch

